Laboratory Shaker (Horizontal Shaking)

MK161

Rotary, elliptical and reciprocate motion

- Compact, space saving design
- Changeable rotary, elliptical and reciprocate motion for mixing, extracting and stirring of samples
- Stable and high torque shaking power and speed with the DC brushless motor
- Shaking frequency and timer are dial setting and digital display
- Shake pause function, timer function and constant operation by one switch
- Selectable mixing, extracting and stirring patterns when used with different shaking stage and racks (optional item)
- Can be placed inside IN602CSW incubator for shaking incubation

Specifications

Model	MK161
Shaking mode	Rotary, Elliptical and reciprocate (Manual operation)
Shaking range	Rotary:30mm Reciprocate: 30mm
Shaking frequency	20~200rpm
Frequency controller	Dial setting, Digital display
Timer	Dial setting, Digital display / Digital 0.1min. (6sec.) to 99.9hr.
Shaking stage dimensions	Main unit : W300 x D254mm, Stage : W290 x D250mm
External dimensions	W350 x D300 x H150mm
Weight	~15kg
Power source 50/60Hz	AC115V 0.5A / AC220V 0.3A



Example of using mounting stage and erlenmeyer flask holder clamps (optional)



Incubator IN602CW slide shaker stage (optional) and MK161 installation example

*Glassware not included.

Operational accessories

Mounting stage



Capacity	Number of erlenmeyer flask clamp
100ml	10pcs
200ml	9pcs
300ml	5pcs
500ml	4pcs
1,000ml	2pcs
Product code	232061

Erlenmeyer flask holder clamp



Product code	Capacity	No. of clamps
232063	200ml	9 pcs.
232062	100ml	10 pcs.
232064	300ml	5 pcs.
232065	500ml	4 pcs.
232066	1,000ml	2 pcs.

^{*}Mounting stage sold seperately

Diagonal rack holder



	Diagonal erlenmeyer flask holder		
Product code Capacity N		No. of unit	
	232067	100ml	3 pcs
	232068	200ml	2 pcs
	232069	300ml	2 pcs

*Mounting	stage	sold	seperately

Diagonal test tube holder			
Product code	Diameter	No. of unit	
232080	ø12mm	50 pcs	
232081 ø16.5 20 pcs		20 pcs	
232082 ø18 20 pcs			

*Mounting stage sold seperately

Diagonal centrifugal tube holder



For spitz tube

1 of Spitz tube			
	Prouct code	Size	No. of units
	232070	15ml	12 pcs.

For 50ml centrifugal tube

For John Centinugal tube		
Product code	Diameter	
232083	ø29mm	

^{*}Mounting stage sold seperately

Non-skid sheet



Product code	Dimension (W×D×H)
232084	290×250×30

^{*}Mounting stage sold seperately

Single spring shaking rack



Dimension (\	N×D×H)		
290×250×66	mm		
Number of test tube:			
ø16mm test	ø16mm test tube×64 (45°inclination)		
Number of e	rlenmeyer flask		
50ml×20pcs, 100ml×10pcs, 200ml×9pcs,			
300ml×5pcs, 500ml×4pcs. 1000ml×2pcs			
Product code	232050		
*Mounting sta	ge not necessary		

^{*}Mounting stage not necessary
This can be set directly to the main
unit.

Two layer spring shaking rack



Dimension (W×D×H)		
290×250×110mm		
Number of test tube:		
ø16mm test tube×64 (45°inclination)		
Number of erlenmeyer flask		
50ml×20pcs, 100ml×10pcs, 200ml×9pcs,		
300ml×5pcs, 500ml×4pcs. 1000ml×2pcs		
Product code	232056	

*Mounting stage not necessary
This can be set directly to the main
unit.

^{*}Glassware not included.